

ePascal Quick Start Guide



- | | | | |
|------------------------|-------------------------------------|---------------------------------------|--------------------|
| 1. Power Button | 2. Emergency Power Supply | 3. Search/Review | 4. Patient Info |
| 5. Probe & Mode Toggle | 6. Custom Key | 7. Cursor | 8. Adjustment Knob |
| 9. M-Mode | 10. Tissue Stiffness Quantification | 11. D-Mode | 12. C-Mode |
| 13. B-Mode | 14. Touch Bar | 15. Body Marker Annotation | 16. Update(Toggle) |
| 17. Measure | 18. Confirm Button (Left) | 19. Confirm Button (Right) | 20. Video Capture |
| 21. Trackball | 22. Image Save | 23. Volume Control | 24. Zoom/Angle |
| 25. Focal Position | 26. Depth Adjustment | 27. One-Touch Optimization (Reserved) | |
| 28. Freeze | | | |

Start

New Patient Information & Probe Selection :

Click the [Patient] button to enter patient information

Click the [Probe] button to select the probe and examination mode

Image Acquisition

B-Mode: Click the [B] button to enter B-mode.

C-Mode: Click the [C] button to enter Color mode. Use the touchscreen and confirmation buttons to adjust the position and size of the sample box.

D-Mode: Click the [D] button to enter PW mode. Use the touchscreen and confirmation buttons to adjust the position and size of the Q-box. Click the [Update] or [D] button to acquire the spectrum.

E-Mode: Click the [E] button to enter Shear Wave Elastography mode. Use the touchscreen and confirmation buttons to adjust the position and size of the Q-box. Click the [Set] button to obtain parameter measurements.

M-Mode: Click the [M] button to enter M-mode. Use the trackball to adjust the sampling line position.

Measurement

Click the [Measure] button to enter application measurement.

Select the measurement tool on the touchscreen, then perform measurements on the target area.

Use the [Erase] button on the touchscreen to clear measurement results from the screen.

Body Marker Annotation

Click the [ABC] button to:

Select the body marker position and confirm probe orientation by rotating the [Angle] knob.

Choose the "Annotation" tab on the touchscreen, then select target location and

Image/Video Storage

add annotations.

Image Storage: Click the [Image] button to save a single frame. Click the [Freeze] button to review and save the target frame.

Video Storage: Click the [Freeze] button to pause imaging, then save dynamic video (forward storage).

Click the [Clip] button for direct video storage (backward storage).